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096083

TA-163.041-KNX02

Instructions for Use for the ZU6985 4-Electrode Conductivity Sensor

The following safety instructions contain the necessary information for the safe use of the product. For any questions relating to safety, simply contact Knick Elektronische Messgeräte GmbH & Co. KG.

1 Safety

Hazards due to temperature or aggressive media are possible, depending on the location of use.

1.1 Intended Use

The sensor is used to measure conductivity and temperature in liquid, aqueous media. The sensor was designed for discontinuous use with Model 703 laboratory conductivity meter and with portable analyzers (COND) in the Knick Portamess and Portavo series.

It is highly resistant to chemicals (materials: glass, platinum, PTFE, FKM). The tube and O-rings can be replaced without changing the cell constant. The sensor is not suitable for sterilization/autoclaving.

1.2 Personnel Requirements

Personnel must receive instruction in how to handle this product.

2 Commissioning

- On unpacking, check the sensor for mechanical damage.
- Connect the sensor cable to the Model 703 conductivity meter (terminal for 4-electrode sensor) or use a ZU0290 adapter (not included) to connect it to a Portamess COND or Portavo COND analyzer.

3 Operation

3.1 Measurement

To obtain a correct measurement, immerse the sensor in the process medium until the liquid level is between the two markings – see dimension drawing (e). Otherwise incorrect measurements can occur. Moving the sensor up and down repeatedly before measurement ensures stable measurement results due to thorough mixing and temperature equalization at the electrodes.

3.2 Temperature Detector

The integrated temperature detector is intended for measuring the temperature of the medium and for conductivity temperature compensation in the analyzer.

4 Maintenance and Cleaning

The sensor can be cleaned using commonly used chemical cleaning agents (acids, alkalis, solvents) and mechanical aids such as brushes or sponges. The tube can be carefully pulled off while rotating it slightly.

CAUTION! Electrode carrier may break. When pulling off or pushing on the tube, do not tilt it. Do not apply mechanical loads to exposed electrode carriers.

5 Calibration

You should calibrate the sensor to obtain precise measurements and for checking the function. The following calibration standards can be used:

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Refer to the user manual of the measuring device for information on how to calibrate.


Specifications

Measuring range	1 µS/cm ... 1 S/cm
Cell constant measured in a saturated NaCl solution	Approx. 1.15/cm, for actual value see nameplate (i)
Temperature	-20 ... 100 °C
Temp detector	Pt1000
Immersion depth	Min. 60 mm Max. 80 mm See dimension drawing (e)
Materials, Dimensions	See dimension drawing legend
Connection	DIN plug, 8-pin (k)

Pin Assignment

1	(+) Current
2	GND
3	(-) Current
4	(+) Voltage
5	(-) Voltage
6	Pt1000
7	Pt1000
8	n.c.

6 Dimension Drawing

Dimension Drawing Legend

a	Cable: length 1.2 m
b	Body: PTFE
c	O-ring: FKM
d	Vent
e	Markings: immersion depth
f	Tube: calibrated precision glass
g	Electrode carrier: glass with platinum electrodes
h	Temperature detector
i	Nameplate: cell constant c, serial number
k	DIN plug, 45325, 8-pin

7 Disposal

The local codes and regulations must be observed when disposing of the ZU6985 4-electrode conductivity sensor.

Knick >

User
Manual

4-Electrode
Conductivity Sensor
ZU6985

Read before commissioning.
Keep for future use.

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